

## Hill-Rom®

### Tri-Flex™ II Bariatric Bed



#### At a glance

- Multiple positions:
  - Single touch Trendelenburg/Reverse Trendelenburg functions
  - Cardiac chair position
- Expandable width frame
- Drop down foot rail
- Independent head and knee positioning
- Scale with easy to read display
- “Fold and roll” storage capability
- Foam mattress with 2 interchangeable 5.5” side bolsters
- Optional Low Air Loss
- Optional trapeze

#### Dependable and Economical

The Tri-Flex™ II Bariatric Bed is a dependable, economical bed that provides support for both patient and caregiver. Features such as expandable width frame, drop down foot rails, battery back up, CPR mechanism and choice of surfaces makes this a reliable, basic choice for the care of your bariatric patients.



**Hill-Rom**

Enhancing Outcomes for Patients and Their Caregivers.™

## Hill-Rom®

### Tri-Flex™ II Bariatric Bed

Tri-Flex bed is available with Synergy® Air Elite surface.

#### Synergy® Air Elite

- Two therapy choices:
  - Low air loss.
  - Alternating pressure.
- Flexibility to customize therapy throughout the healing period.
- Individual air cushions:
  - Designed for maximum pressure relief.
  - Connected to side panel with CPC connector for easy detachment and cleaning.
- Polyurethane fitted cover reduces shear and friction.
- Multiple function choices powered by blower technology allow mattress to perform rhythmic inflation and deflation to redistribute and reduce patient interface pressure.
- 8 digitally-selected scales of soft/firm comfort control.
- Vented air cushions provide effective low air loss therapy.
- Waterproof and vapor permeable ticking.



- Auto firm function provides uniform firmness to accommodate caregiver procedures.
- Safety features:
  - Power failure triggers auto alarm.
  - Main power switch can be shut off in an emergency.
  - CPR hose disconnects for rapid deflation.
  - 2.5" foam base provides support in the event of a power failure.
- Double insulated blower provides quiet operation.
- 39" - 500 lbs., 48" - 700 lbs.



Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2008 Hill-Rom Services, Inc ALL RIGHTS RESERVED.

145598 rev2 2/28/08

## Technical Specifications

### Tri-Flex II

Low position - sleep deck to floor.....	14" (36 cm)
High position - sleep deck to floor.....	28" (71cm)
Maximum overall width.....	40.75" (109 cm)
Maximum overall length.....	93.875" (939 cm)
Sleep deck width.....	37 - 48" (94 - 122 cm)
Sleep deck length.....	86" (218 cm)
Fold and roll storage (H x W x L).....	74.5" x 32" x 46.5" (89 x 81 x 118 cm)
Casters.....	6" (15 cm) steer and lock
Weight Capacity.....	1,000 lbs. (453.6 kg)
Safety listing.....	UL60601-1
Options.....	Scale (Standard with rental) Trapeze Low airloss surface

# Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

Capital **800-445-3730**  
Rental **800-638-2546**  
Canada **800-267-2337**

[www.hill-rom.com](http://www.hill-rom.com)