

Hill-Rom®

TotalCare® Certified Reconditioned Beds



Quality. Safety. Value.



No one knows Hill-Rom beds better than Hill-Rom does! With our state-of-the-art reconditioning process, all TotalCare® beds meet UL standards and are assembled in an ISO certified process.

Plus, enjoy many of the same features of a new TotalCare® bed on a certified reconditioned bed!

Hill-Rom®

Enhancing Outcomes for Patients and Their Caregivers.™

TotalCare® Bed Reconditioning Process

Each bed goes through a series of inspections, predetermined repairs and component replacements.

Beds are thoroughly hand cleaned and painted using matching color touch-up paint.

All frames are upgraded to meet HBSW requirements.

All internal components are replaced with a "non-mercury" equivalent (as needed).

Each bed is inspected by Hill-Rom's quality inspectors using a computerized final inspection process.

170 safety and quality checks are performed, including the same functionality checks that a new bed would receive!

We replace:

- Surface
- All batteries
- Accessory outlet
- Footboard
- Night light bulb
- Hydraulic power unit (if needed)
- All high wear touch points and worn areas
- Casters
- Blower
- Headboard
- Missing or damaged labels
- Power cords
- Graphical Caregiver Interface (GCI®)
- Software

We Upgrade

- All beds equipped with a GCI® (graphical caregiver interface) will include a 30-degree "head of bed" (HOB) alarm.
- Beds equipped with a "nurse call" option will include universal TV (UTV) as well as the NaviCare® Patient Safety Module compatible communication PC board.
- Beds equipped with a scale are upgraded to include high capacity load beams and corresponding scale PC board.
- TotalCare SpO₂RT® beds with low air loss, will have a 550 lb safe working load.
- Siderails to meet HBSW recommendations.

TotalCare® bed options & accessories

Any of the great options available on a new TotalCare® bed are also available on a Certified Reconditioned TotalCare bed! Choose from options such as:

- Nurse call
- NaviCare®
- Line managers
- IntelliDrive®
- Low air loss
- And more!

Hill-Rom's turning green!

At Hill-Rom, we do our part to be good stewards of the environment by:

- Removing mercury from Hill-Rom beds and competitor beds
- Recycling plastic
- Properly disposing of circuit boards and batteries
- Using powder and water based coatings
- Recycling cardboards, plastics, metals and office paper
- Recycling and reclaiming glass, metals and plastics for used fluorescent lighting, batteries and circuit boards

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2008 Hill-Rom Services, Inc ALL RIGHTS RESERVED.

154842 rev 1 9/19/08

Hill-Rom®

Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com