

Hill-Rom®

Non-Powered Surface Portfolio



		NP50	NP100	Tempur-Pedic®	AccuMax Quantum™ VPC	AccuMax Quantum™ Convertible	Hill-Rom® NP200	NP50 Bariatric	NP100 Bariatric
Surface L x W x D (in)		80 x 35.5 x 6	80 x 35.5 x 6	80 x 35.5 x 6	80 x 35 x 7	80 x 35 x 7	80 x 35.5 x 6	80 x 39 x 7	84 x 38 x 7
Additional Lengths		76 and 84	76 and 84	84	75, 84 and 86	84, 86	84, 86	84	N/A
Additional Widths			32, 35 and 38.5	35 and 36	32	32	N/A	48	38.5
Features	Benefits								
Frame Compatibility:									
- Flat Deck		X	X	X	X	X	X	X	X
- VersaCare®			X	X	X	X	X		
- TotalCare®			X	X	X	X			X
- CareAssist®			X	X	X	X	X		
- Stretchers				X	X				
Psychiatric Surfaces		X							
Shear liner	Reduce shear forces which cause pressure ulcers		X	X			X		X
High Quality Ticking Material	Reduces friction which can lead to pressure ulcer formation, provide stretch to enhance immersion and redistribute pressure	Nylon	Bidirectional stretch	Bidirectional stretch	Nylon*	4-way stretch top ticking (Optional Nylon is available)	Bidirectional stretch	Nylon	Bidirectional stretch
Enhanced Heel Zone	Reduces and redistributes pressure		X	X	X	X	X		X
Multi-zone construction	Surface designed and constructed to match the body's proportional differences thereby maximizing redistribution		X	Muti-zoned, Multi-layered	Muti-zoned, Multi-layered	Muti-zoned, Multi-layered	X		X
Fluid Impervious RF Welded Seams	Surface is constructed to thwart cross contamination between patients		X	X			X		X
Latex/DEHP free	Reduces possible allergic reactions	X	X	X	X	X	X	X	X
PVC free	Reduces possible allergic reactions		X	X			X		X
Non-flammable	Meets Fire Safety Standards Cal 129, CPSC 16CFR1633 and Boston IX-11	X	X	X	X	X	X	X	X
Designed to reduce entrapment when used with recommended frames	Helps provide a safe environment for patients	X	X	X	X	X	X	X	X
Firm Perimeter	Provides stable surface for egress	Bolstered Versions Available	X	X	X	X	X	X	X
Optional Powered Therapy						Alternating Pressure Therapy			
Compatible with bed exit alarms	Enhances synergy between systems	X	X	X	X	X	X	X	X

*Optional 4-way stretch available



Enhancing Outcomes for Patients and Their Caregivers.™

Hill-Rom®

Non-Powered Surface Portfolio



		NP50	NP100	Tempur-Pedic®	AccuMax Quantum™ VPC	AccuMax Quantum™ Convertible	Hill-Rom® NP200	NP50 Bariatric	NP100 Bariatric
Features	Benefits								
Mattress (foam components) compressibility tested life	Enhances product life	5 years	5 years	10 years	5 years	5 years	5 years	5 years	5 years
Warranty	Provides for reimbursement	2 years	5 years	10 years**	5 years limited; 3 years prorated on Air System; 2 years non prorated on top and bottom covers and Foam Containment System	5 years limited; 3 years prorated on Air System; 2 years non prorated on top and bottom covers and Foam Containment System; 1 year on powered pump option	5 years	2 years	5 years
Weight limit (lbs)	Able to accommodate greater number of patients	350	350-Standard 500-VersaCare®	500	500	500	500	700	500
Mattress weight	Allows patient handler to easily maneuver surface	24 lbs	21 lbs	35 lbs	29-34 lbs	29-34 lbs	20 lbs	35 lbs	26-30 lbs
Handles	Allows patient handler to easily maneuver surface		Flat Deck, CareAssist®	Flat Deck, CareAssist®	X	X	Flat Deck, CareAssist®		X
Wipe down surface material	Easy cleaning	X	X	X	X	X	X	X	X
Anti-microbial foam	Controls bacterial growth	X	X	X	X	X	X	X	X
Anti-microbial Ionized Silver Ticking option*	To inhibit growth of bacteria that may affect this product		X	X			X		X
Anti-static, acid resistant and waterproof	Extends product life	X	X	X			X	X	X
Mattress retention system or anti-skid bottom	Helps keeps a mattress in place on the bed	X	X	X			X	X	X
Heavy duty mattress material	Enhances mattress durability	X			X	X		X	
Available replacement covers, air cells, components	Extends product life		X	X	X	X	X		X
Clinical positioning		Prevention	Prevention	Prevention and treatment, Stage I, II	Prevention and treatment, Stage I and II and uncomplicated III's and IV's with multiple turning surfaces	Prevention and treatment, Stage I and II and uncomplicated III's and IV's with multiple turning surfaces	Prevention and treatment, Stage I and II and uncomplicated III's and IV's with multiple turning surfaces	Prevention	Prevention

*Use of this product does not reduce or alter hospital infection protocol. **Requires replacement ticking from Hill-Rom within 5 years.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2009 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

154936 rev 2 12/01/09



Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com