

Hill-Rom®

Basis™ Headwall



Tired of placing institutional headwalls in your hospitality-inspired patient rooms?

The Basis™ headwall system was designed to help you quickly and easily integrate a headwall, with all of its mission-critical services, into a hospitality-inspired patient room design.

A clean, linear design free of faceplates and painted-metal surfaces helps you create an environment that evokes hospitality while still delivering an environment that supports superior patient care.

More than just a pretty headwall, the Basis headwall supports your hospital's lean initiatives.

The Basis headwall's superior adaptability makes it an obvious choice for facilities adopting lean processes in the patient care environment. Traditional headwalls use copper tubing to deliver gas service. The Basis headwall uses flexible hoses, like those used in the operating room environment for years, so you can add or relocate a service to optimize a workflow within minutes, instead of hours or days. And, you don't have to take the room out of service!

Hill-Rom®

Enhancing Outcomes for Patients and Their Caregivers.™

Hill-Rom®

Basis™ Headwall



Helping you address HIPAA/privacy concerns

Patients and their families deserve an environment where they can openly and candidly communicate, free from the fear of being overheard by (or overhearing) their “neighbors.” A headwall is an integral piece of the wall design, and the Basis headwall is the only headwall designed and tested to meet and exceed the American Institute of Architect’s recommended sound transmission class (STC) ratings.

Sustainable design... caring for the patient and the community

The Basis™ system is the only headwall system that ships complete with a “cradle-to-cradle” lifecycle strategy. The “cradle-to-cradle” design philosophy embraces the idea that at the end of a product’s useful life, it can be turned into something else to close the cycle so that ultimately there is no waste. The Basis headwall helps you care for your patient and your community.

For more information, contact your Hill-Rom representative at 800-445-3730 or visit www.hill-rom.com/basis.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2010 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.



Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com