


Hill-Rom Power Column, Flatwall Headwall System and Contour Headwall System



Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
© 2007 Hill-Rom Services Inc.
ALL RIGHTS RESERVED.

CTG019ra rev 2 1/8/07

USA • 800-445-3730 • Fax 812-934-8189
Canada • 800-267-2337 • Fax 905-206-0561
www.hill-rom.com


Hill-Rom



Enhancing Outcomes for Patients and Their Caregivers.™

Hill-Rom




Enhancing Outcomes for Patients and Their Caregivers.™







	Power Column	Flatwall	Contour	Wall Mount
 <p>P976-02 Infusion Support Arm</p> <ul style="list-style-type: none"> Designed to hold two transfer poles (P158A) and can be located almost anywhere on the track 	●	●		
 <p>P984-00 Infusion Pump Holder</p> <ul style="list-style-type: none"> Four hangers, requires two P970A01 universal holders 	●	●		
 <p>P972-01 Post-Type Utility Mount</p> <ul style="list-style-type: none"> Requires P970A01 Universal Holder Length of post - 2 1/4" 	●	●		
 <p>P973-01 Utility/Bottle Slide Track Mount</p> <ul style="list-style-type: none"> Requires P970A01 Universal Holder 	●	●	●	
 <p>P1020D1 Utility/Bottle Slide Flush Mount</p>	●	●	●	
 <p>P971H00 Slide-Mounted IV Hanger</p> <ul style="list-style-type: none"> Outer Dimensions: 2"W x 4 1/2"H Inner Dimensions: 1 3/4"W x 4"H Requires P973-01 or P1020D1 Utility/Bottle Slide Track Mount Requires P970A01 Universal Holder 	●	●	●	








	Power Column	Flatwall	Contour	Wall Mount
 <p>P985A01 23" Power Column Accessory Bar</p> <p>For accessories that mount on this arm, please see Integris B accessories brochure</p>	●			
<p>P985A05 15" Power Column Accessory Bar</p> <p>For accessories that mount on this arm, please see Integris B accessories brochure</p>	●			
<p>P985A02 19" Flatwall Accessory Bar</p> <p>For accessories that mount on this arm, please see Integris B accessories brochure</p>		●		
<p>P985A03 28" Flatwall Accessory Bar</p> <p>For accessories that mount on this arm, please see Integris B accessories brochure</p>		●		
<p>P985A04 37" Flatwall Accessory Bar</p> <p>For accessories that mount on this arm, please see Integris B accessories brochure</p>		●		
<p>P985A06 Contour Accessory Bar</p> <p>For accessories that mount on this arm, please see Integris B accessories brochure</p>			●	

	Power Column	Flatwall	Contour	Wall Mount
 <p>P979-00 6' Vertical Wall Tracks with Keyway at Center Wall mount only</p>				●
 <p>P985-15 15 1/8" Power Column Accessory Bar For accessories that mount on this arm, please see Horizon/Uniflex accessories brochure</p>	●			
<p>P985-11 23" Power Column Accessory Bar For accessories that mount on this arm, please see Horizon/Uniflex accessories brochure</p>	●			
<p>P985-12 19" Flatwall Accessory Bar For accessories that mount on this arm, please see Horizon/Uniflex accessories brochure</p>		●		
<p>P985-13 28" Flatwall Accessory Bar For accessories that mount on this arm, please see Horizon/Uniflex accessories brochure</p>		●		
<p>P985-14 37" Flatwall Accessory Bar For accessories that mount on this arm, please see Horizon/Uniflex accessories brochure</p>		●		
<p>P985-16 Contour Accessory Bar For accessories that mount on this arm, please see Horizon/Uniflex accessories brochure</p>			●	



	Power Column	Flatwall	Contour	Wall Mount
 <p>P158A Transfer Pole <ul style="list-style-type: none"> • Includes four IV hangers • Retracted Position - 32" • Fully Extended - 57 1/4" • Diameter of top pole - 1/2" • Diameter of outer pole - 1" • Requires P976-02 Infusion Support Arm </p>	●	●		
 <p>P972L01 Post-Type Long Utility Mount <ul style="list-style-type: none"> • Requires P970A01 Universal Holder • Length of post - 12" </p>	●	●		
 <p>P970A01 Universal Holder </p>	●	●	●	
 <p>P985P03 Accessory Arm (left-shown) P985P04 Accessory Arm (right) <ul style="list-style-type: none"> • Right Hand 18" • Left Hand 18" For accessories that mount on these arms, please see Horizon/Uniflex Accessories brochure</p>	●	●	●	●






	Power Column	Flatwall	Contour	Wall Mount
<p>P962CMU Hewlett-Packard Merlin CPU/Satellite Rack Mount</p> <ul style="list-style-type: none"> • Optional mounting of CPU on back of Power Column 	●			
 <p>P962-00 Tilt Swivel Monitor Arm</p> <ul style="list-style-type: none"> • Arm requires built in channel support brackets inside of headwall • Requires monitor mfg. name and model # 	●	●		
 <p>P997L00 Light Weight Monitor Arm</p> <ul style="list-style-type: none"> • Arm requires built in channel support brackets inside of headwall • Requires monitor mfg. name and model # 	●	●		
 <p>P100702 Monitor Slide, Wall Mount Only</p>				●

	Power Column	Flatwall	Contour	Wall Mount
 <p>P967B00 STAT Clock•Timer Panel Mount</p> <ul style="list-style-type: none"> • For replacement only 	●	●		
 <p>P964B00 STAT Clock•Timer Pod Mount</p> <ul style="list-style-type: none"> • Hard wiring required • Requires P970A01 Universal Holder 	●	●		
 <p>P967B01 STAT Clock•Timer Wall Mount</p> <ul style="list-style-type: none"> • Requires four gang masonry box 3 1/2" deep provided by others 	●	●		
 <p>P973-C1 Cable Organizer</p> <ul style="list-style-type: none"> • Requires P970A01 Universal Holder 	●	●		
 <p>P958-00 Floor Mounted Docking Shoe</p>				
 <p>P1020C1 Extended Bottle Slide</p> <ul style="list-style-type: none"> • Dimensions - 1 5/8"W x 4 5/16"H • Requires accessory bar or arm to mount 	●	●	●	

	Power Column	Flatwall	Contour	Wall Mount
 <p>P963C01 Arm Lamp - Incandescent Light</p> <ul style="list-style-type: none"> • Arm extends to 44" • Requires power outlet and P100803 Mount for track mounting • 13" Cord Length 	●	●	●	●
 <p>P963C02 Arm Lamp - Fluorescent and Incandescent Bulbs</p> <ul style="list-style-type: none"> • Arm extends to 44" • Requires power outlet and P100803 Mount for track mounting • 13" Cord Length 	●	●	●	●
 <p>P963M00 - Wall Mount for Arm Light</p>				●
 <p>P100803 - Mount Bracket</p>	●	●	●	
 <p>P966A00 Utility Light</p> <ul style="list-style-type: none"> • Requires power outlet • 24" Cord Length 	●	●	●	
 <p>P1022PH Resuscitator Bag Holder</p> <ul style="list-style-type: none"> • Requires P970A01 Universal Holder 	●	●	●	
 <p>P973-B1 Bird Blender Mount</p> <ul style="list-style-type: none"> • Requires P970A01 Universal Holder 	●	●	●	

	Power Column	Flatwall	Contour	Wall Mount
 <p>P978-01 Aneroid Sphygmomanometer with coil and a standard size cuff</p> <ul style="list-style-type: none"> • The Tycos aneroid sphygmomanometer mount, cuff basket and universal holder are ordered separately. 	●	●	●	
 <p>P975-01 Aneroid Mount and Cuff Basket</p> <ul style="list-style-type: none"> • Requires P970A01 Universal Holder 	●	●	●	
 <p>P975A01 Aneroid Mount for a Welch Allyn Model 7670-01 .</p> <ul style="list-style-type: none"> • Requires a P970A01 Universal Holder 	●	●	●	
 <p>P973-S2 Integrated Diagnostic System Mount/ Welch Allyn Model 767</p> <ul style="list-style-type: none"> • Requires P970A01 Universal Holder 	●	●	●	
 <p>P973-S1 Integrated Diagnostic System Mount/ Welch Allyn Model 747</p> <ul style="list-style-type: none"> • Requires P970A01 Universal Holder 	●	●	●	

	Power Column	Flatwall	Contour	Wall Mount
 <p>P951W03 Large Pivoting Storage Basket</p> <ul style="list-style-type: none"> • Inner dimensions - 9 9/16" W x 9 9/16" H x 12 7/8" D • Functional loading up to 24 lbs. • Light neutral color 	●	●	●	
 <p>P951S03 Small Pivoting Storage Basket</p> <ul style="list-style-type: none"> • Inner dimensions - 9 9/16" W x 9 9/16" H x 4" D • Functional loading up to 4 lbs. • Light neutral color 	●	●	●	
 <p>P952S01 Small Storage Basket</p> <ul style="list-style-type: none"> • Inner dimensions 21 3/8" W x 4 3/8" H x 2 7/8" D • Functional loading up to 3 lbs • Light Neutral Color 	●			
 <p>P952L01 Large Storage Basket</p> <ul style="list-style-type: none"> • Inner dimensions 20 7/16" W x 8 13/16" H x 6 7/8" D • Functional loading up to 15 lbs • Light Neutral color basket with Black bumper 	●			
 <p>P672E03 Storage Basket for 19" Flatwall</p> <ul style="list-style-type: none"> • Inner dimensions: 18 11/16" W x 7" H x 9 1/16" D • Functional loading up to 15 lbs 		●		
 <p>P672E04 Storage Basket for 28" Flatwall</p> <ul style="list-style-type: none"> • Inner dimensions: 27 11/16" W x 7" H x 9 1/16" D • Functional loading up to 20 lbs • Light Neutral Color with Black Bumper 		●		

	Power Column	Flatwall	Contour	Wall Mount
 <p>P986S01BMO7 Small Swivel Shelf</p> <ul style="list-style-type: none"> • Light neutral laminate • Overall dimensions - 15" x 12" • Weight capacity - 50 lbs. 	●	●	●	
 <p>P986-03BMO7 Large Swivel Shelf</p> <ul style="list-style-type: none"> • Light neutral laminate • This swivel shelf can be rotated 360°. • Overall dimensions - 15" x 22" • Weight capacity - 50 lbs. 	●	●	●	
 <p>P961E02BMO7 Adjustable Shelf for Power Column (shown)</p> <ul style="list-style-type: none"> • Light neutral laminate • Overall dimensions - 22 1/8" W x 12 1/8" D • Weight capacity - 50 lbs. 	●*			
<p>P673E21BMO7 Adjustable Shelf for 19" Flatwall</p> <ul style="list-style-type: none"> • Light neutral laminate • Overall dimensions - 22 3/8" W x 12 1/8" D • Weight capacity - 50 lbs. 		●		
<p>P673E30BMO7 Adjustable Shelf for 28" Flatwall</p> <ul style="list-style-type: none"> • Light neutral laminate • Overall dimensions - 29 3/8" W x 12 1/8" D • Weight capacity - 50 lbs. 		●		
<p>P673E16B Adjustable Shelf for Contour</p>			●	
 <p>P674E01 19" Flatwall Guard/Bumper (shown)</p> <p>Also Available:</p> <ul style="list-style-type: none"> • P674E02 28" Flatwall Guard/Bumper • P674E03 37" Flatwall Guard/Bumper • P953-00 23 1/8" Power Column Guard/Bumper • P953-01 15 1/8" Power Column Guard/Bumper • P674E07 Contour Guard/Bumper 		●		
 <p>P674E05 28" Flatwall Bumper with Three Bottle Racks (shown) ...</p> <p>Also Available:</p> <ul style="list-style-type: none"> • P674E04 19" Flatwall Bumper with Two Bottle Racks • P674E06 37" Flatwall Bumper with Four Bottle Racks 		●		

*Not available for the 15 1/8" Power Column