

Hill-Rom®

Power Columns



Power Column

Maximum access, maximum flexibility

At a glance

- Four-sided access to the patient
- Efficient use of premium space
- Two-sided accessory tracks on all four corners of the column
- Numerous accessories available
- Two sizes of power columns (23" and 15") for maximum flexibility

Hill-Rom power columns offer the best of both worlds: 360° access to critically ill patients, and the best use of space in the crowded acute care environment.

Hill-Rom power columns allow caregivers to move freely around the bed, including the critical area around the head of the patient. Free of the wall, the bed can be pivoted around to allow the patient to see out the window or watch television, without disturbing vital services.


Hill-Rom®

Enhancing Outcomes for Patients and Their Caregivers.™

Hill-Rom®

Power Columns

Built and configured just for you



Each power column is designed and built to individual specifications. The unique accessory track system provides full-length, two-sided tracks on all four corners, allowing staff to position equipment where it's needed most.

Both power columns offer full service capacity, making it easy for staff to focus on patient care.

Technical Specifications

Column cross-section

Feature	Dimension
Power Column.....	23.125" W x 7" D
Small Power Column	15.125" W x 7" D

Floor Space (Base)

Feature	Dimension
Power Column.....	21.5" W x 6" D
Small Power Column	13.5" W x 6" D

Ceiling Space

Feature	Dimension
Power Column.....	23.25" W x 10.5" D
Small Power Column	15.25" W x 10.5" D

Height

Floor to Ceiling.....	.8'0" – 9'6" in 1" increments
-----------------------	-------------------------------



Small Power Column
Side View



Small Power Column
Front View

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2008 Hill-Rom Services, Inc ALL RIGHTS RESERVED.

CTG267rb rev 2 10/15/08

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-638-2546
Canada 800-267-2337

www.hill-rom.com