


Hill-Rom Vollman Proner



At a glance

- Facilitates safe movement of a patient to the prone position
- Allows access to all sides of the patient
- Easy to clean
- Supports patients up to 300 pounds (136 kg)

Maintain access while proning immobile patients

With a single, cost-effective tool, caregivers can place patients in the prone position and still maintain access to all sides of the body.

Positioning immobile patients is essential to improving outcomes and avoiding adverse events. The prone position has been proven to help facilitate gas exchange, improve pulmonary blood flow, and increase overall lung volumes.

Vollman Proner provides extra support and stability while three caregivers turn the patient, and then supports the patient in the prone position, leaving the abdomen freely suspended. Procedures such as suctioning are easier to perform.

Vollman Proner allows patients to stay in the same bed as they recover, eliminating a disruptive transfer. Patients can move more smoothly through the progressive mobility process by staying in place.



Therapeutic appearance

Unlike some proners, the Vollman Proner looks gentler and less industrial. Families find its softer appearance easier to accept.




Cost effective

A tool that's priced to own, saving money on rentals. In fact, the Vollman Proner can be purchased for the price of a week's rental of a proning bed.



Easy to clean

Cushions and support rods wipe clean easily for multi-patient use and effective infection control.



A stationary chest support cushion and adjustable pelvic cushion supports the patient. The abdomen is freely suspended for diaphragmatic movement and improved functional residual capacity (FRC).

The prone position helps facilitate gas exchange, improve pulmonary blood flow, and increase overall lung volumes.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
© 2007 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

CTG297 Rev. 2 4/17/07