

**Hill-Rom**  
**UniFlex® 1000**  
**Headwall System**



**At a glance**

- Delivered pre-plumbed and pre-wired to permit single point connections for each service type
- Available with an optional vertical chase
- Ancillary equipment is easily attached along the upper and lower accessory tracks
- Comes in variable lengths and a wide choice of laminate colors

**A practical solution to equipment and service management**

The UniFlex® 1000 Headwall (Part #P1002), from Hill-Rom, brings a clean, contemporary look to intermediate care areas of the hospital.

This efficient headwall provides an alternative to in-wall construction systems. It is delivered pre-plumbed and pre-wired to permit single point connections for each service type and mounts directly to the wall with brackets that are provided.

**Hill-Rom®**

Enhancing Outcomes for Patients and Their Caregivers.™

# Hill-Rom UniFlex® 1000 Headwall System



Available with an optional vertical chase, the UniFlex® headwall organizes vacuum, oxygen, air, electrical and communications within easy reach of the staff. Ancillary equipment is easily attached along the upper and lower accessory tracks allowing the system to adapt to changing care requirements.

The headwall comes in variable lengths and a wide choice of laminate colors. It may be used alone, or as a system with the addition of lighting and a Bed Locator from Hill-Rom.



Integris®  
Patient Light  
(P6961L)



Bed Locator

## ▶ Technical data

<b>Length (1 rail)</b> .....	3' minimum up to 9' maximum
.....	(91.44 cm up to 274.32 cm)
.....	For units longer than 9', rails easily connect together
<b>Rail Height</b> .....	14" (35.56 cm)
<b>Rail Depth</b> .....	3-15/16" (10 cm)
<b>Rail Distance from floor to the device centerline</b>	
standard .....	48" to 60" (121.92 cm to 152.4 cm)
<b>Vertical Chase</b> .....	5-5/8" deep by 13-3/4" wide
.....	(14.29 cm x 34.93 cm)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2006 Hill-Rom Services, Inc.  
ALL RIGHTS RESERVED.

CTG309 Rev. 2 8/15/06

USA • 800-445-3730 • Fax 812-934-8189  
Canada • 800-267-2337 • Fax 905-206-0561  
www.hill-rom.com

**Hill-Rom**

Enhancing Outcomes for Patients and Their Caregivers.™