

Hill-Rom®

Art of Care® 3-Position Recliner



Metropolitan

Opera

Aero/Artisan



Perfect Tailoring: The recliner is tailored with style and perfection in mind.

At a glance

- Exceptionally smooth and easy to recline mechanism
- Intuitive and user friendly operation without cumbersome lever
- Front casters with total lock, positioning handles and arm grips to facilitate user handling
- Accessory hooks on both sides
- STEELock™ construction ensures superior structural integrity
- Available with 21" seat height
- Optional removable arm to facilitate lateral ingress and egress
- Optional accessory holder and tray table

Hill-Rom® 3-Position Recliner delivers style and a relaxing experience for visitors and patients.

Created with user experience in mind, Hill-Rom's 3-Position Recliner integrates with our four Art of Care styles: Artisan™, Metropolitan™, Aero™ and Opera™. Using spring box seat construction for seating comfort and an exceptionally smooth and easy to use mechanism, the recliner construction stands for cleanability and durability. Positioning handles, locking casters, accessory hooks, optional removable arm and tray table are just a few of the distinctive features.

Hill-Rom®

Enhancing Outcomes for Patients and Their Caregivers.™

Hill-Rom®

Art of Care® 3-Position Recliner



Accessory Hooks

Located on both sides, the spring loaded hooks retract automatically in place for cleanliness and safety.



Positioning Handles

On the back of the seat, the ergonomically placed handles allow for easy positioning.



Cleanout Space

The convenient cleanout space makes cleaning simple.



Smooth and Protected Mechanism

User friendly and smooth, the 3-position recline mechanism is built for durability. For increased safety, the integrated ottoman mechanism is protected.



Arm Grip

Extended arm grips assist user during ingress and egress.



Total Lock Casters on Front

The total lock casters ensure stability or smooth movement on demand.



Optional Tray Table

Conveniently placed on any side, the tray table rotates easily out of the way. It can be stored on the side of the recliner and lift back in place in one movement. Additionally, it can be removed to reveal an accessory holder.



Optional Removable Arm

Located on either side or both, the removable arm allows for easy lateral ingress and egress.

Specifications

3-Position Recliner



R3SR	R3SM	R3SO	R33R	R33M	R33O
Aero/Artisan	Metropolitan	Opera	Aero/Artisan	Metropolitan	Opera

Width: 28.5" **Depth:** 39.75" **Height:** 49.5" 3-Position Recliner with 21" seat
SW: 23.5" **SD:** 22" **SH:** 19.25"

COM Ydg. 6.25

Maximum Weight: 300 lbs.

Product	CAL 133	Options	Caps Finish	Upholstery
PLEASE SPECIFY THE FOLLOWING FOR EACH 3-POSITION RECLINER				
R3SR, R3SM, R3SO or R33R, R33M, R33O	CAL 133	Right tray table Left tray table Right Removable Arm Left Removable Arm Both Arms Removable	See Hill-Rom® Finish Options brochure ¹ or visit www.hill-romdesign.com	See Hill-Rom® Performance Collection brochures ² or visit www.hill-romdesign.com

1. CTG 255rb
2. CTG 659

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2008 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

135064 rev 4 7/14/08



Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com