

Hill-Rom®

Art of Care® 4-Position Recliner



1 **As shown**

*R4SP: Fabric; Silverado Ginger (L74)
Finish Honey Maple*



2



3



4

At a glance

- Four Positions – seating (1), TV (2), recline (3) and feet above head (4)
- Exceptionally smooth and easy to recline mechanism
- Front locking casters, positioning handles and arm grips to facilitate egress/ingress
- Accessory hooks on both sides
- SteelLock™ construction ensures superior structural integrity
- Optional removable arm, tray table or accessory bracket on left, right or both
- Choose from either wood or urethane arms

Hill-Rom's 4-Position Recliner delivers style, comfort and functionality for patients and caregivers.

Created to provide an aesthetically pleasing design that will enhance the overall experience for the patient and visitor. Hill-Rom's 4-Position recliner integrates with our four Art of Care® styles. Using spring box construction for seating comfort and an exceptionally smooth and easy-to-use mechanism. This recliner provides performance, durability and utility that you have come to expect from Hill-Rom. Positioning handle, front locking casters, accessory hooks, optional removable arm, tray tables, urethane arms caps and accessory brackets for an IV pole or O₂ tank are just a few of the distinctive options.

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

Hill-Rom®

Art of Care® 4-Position Recliner

Hill-Rom's 4-Position Recliner has standard features or options that will ensure a comfortable, functional and stylish sitting surface.



Accessory Hooks

Located on both sides, the spring loaded hooks retract automatically in place for cleanliness and safety.



Smooth and Protected Mechanism

User friendly and smooth, the 4-position recliner mechanism is built for durability. For increased safety, the integrated ottoman mechanism is protected.



Arm Grip

Extended arm grips assist user during ingress and egress.



Optional Accessory Bracket

For IV pole or O₂ tank.



Optional Tray Table

Conveniently placed on any side the tray table rotates easily out of the way to be stored.



Optional Removable Arm

Located on either side or both, the removable arm allows for ingress and egress.

Product	CAL 133	Options	Caps Finish	Upholstery
PLEASE SPECIFY THE FOLLOWING FOR EACH 4 POSITION RECLINER				
R4SP	Cal 133 available	Right or left tray table. Right, left or both removable arm. Right, left or both accessory bracket. Urethane arms.	See Hill-Rom. finish options brochure or visit www.hill-rom.com	See Hill-Rom. fabric cards or visit www.hill-rom.com

Visit www.hill-rom.com to "Design Your Own" Art of Care® 4-Position Recliner

Technical Specifications



R4SP 4-Position Recliner

W 29.00 D 40.00 H 50.00 SW 23.00 SD 22.00 SH 19.00

Footprint Fully Reclined

W 29.00 D 69.00

COM Ydg: 6.25

Maximum Weight: 300 lbs

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2009 Hill-Rom Services, Inc ALL RIGHTS RESERVED.

151602 rev 2 8/17/09



Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com