

Hill-Rom®
 Advanta™ 2 Bed



Designed for patient safety and caregiver safety/satisfaction.

Top Advantages

Hill-Rom
 Advanta™ 2 Bed

Chair Position: Caregivers can place patient in a chair position with the touch of one button. Bed remains at low height to help reduce risk of patient falls.	Yes
Low Height: Helps reduce the risk and impact of patient falls.	15.75"
Built-in Battery Backup: Backup for all bed articulations.	Yes
Point-of-Care® Controls: All caregiver controls are conveniently located in a central location, on both sides of the bed, where patient and caregiver interaction occurs.	Yes
DC Motors: Provides quiet and smooth bed operation.	Yes
Brake Alarm: Alerts the caregiver when the bed is plugged in and the brake is not set.	Yes
Care Grip: Ergonomically correct handle to grip for safe bed entry/exit reducing fall risk.	Yes
Digital Head of Bed Indicator: Located conveniently in the siderail. Large display, clear indication of head angle.	Yes
Electric Vascular Foot: Eliminates the need for the nursing staff to lift the weight of the bed and patient's legs.	Yes
Integrated Bed Extender: Retractable foot mechanism (surface 84" - 87").	Yes



Hill-Rom®

The Advanta™ 2 Bed: Designed for patient safety and caregiver safety/satisfaction

Hill-Rom Features on STANDARD FRAME (NO ADDITIONAL COST)	Hill-Rom Advanta™ 2 Bed
Low Height: Lowest.	Yes
Dampened Emergency CPR: Dampened CPR gently lowers patient when activating CPR.	Yes
Brake Alarm: Alerts caregiver when brakes are not set.	Yes
Yellow LED Nightlight: Always on; located on both sides at patient egress point.	Yes
Drop-Down Siderails: One handed straight up/down siderail release.	Yes
Non-Moving Headboard: Does not move with hi/lo motion of bed; protects headwalls and lights.	Yes
Built-in Battery Backup: Backup for all bed articulations.	Yes
Point-of-Care® Controls: Siderail controls give caregiver electric adjustment for height, head, knee, foot, Trend/Reverse Trend, chair, scale and bed exit (no controls at foot of bed).	Yes
Patient Pendant Controls: Easy access patient controls in a pendant. Can be stowed in upper and lower siderails.	Yes
Line-of-Site® Angle Indicator: Head of bed and Trendelenburg indicators built into siderails.	Yes
Bed-Not-Down Indicator: Indicator in siderail lit when bed's not in lowest position.	Yes
Electric Vascular Foot: Eliminates the need for the nursing staff to lift the weight of the bed and patient's legs.	Yes
Chair Position: Head of bed angle 65 degrees; maintains low 15.75" height.	Yes
Integrated Bed Extender: Retractable foot mechanism (surface 84" - 87").	Yes
Obstacle Safeguard™ Feature: Floating foot section protects against damage to equipment and furniture.	Yes
Brake/Steer Pedals: Located at all four corners.	Yes
Brake Lock: Locks all four casters and prevents bed from rolling and moving sideways.	Yes
Shearless Pivot® Mechanism: Maintains optimal patient placement while reducing shearing forces when head of bed is raised.	Yes
Ergonomically-Designed Push Handles	Yes
Cord Wrap Clip / IV Pole Storage: Design upgrade for cord management.	Yes
Integrated Line Management: Industry's first line managers integrated into siderail for effective line organization.	Yes
Customizable Graphics	Yes
Bed Identity Option: Customize the bed with color, graphics and unit identification to compliment your facility aesthetics and asset management.	Yes
30 Degree Head of Bed Alarm	Yes
Patient Storage: Industry's first hospital bed with convenient in-rail storage for personal items.	Yes
Care Grip: Ergonomically correct handle to grip for safe bed entry/exit reducing fall risk.	Yes
Four Foley Bag Hooks	Yes
IV Pole Storage: Positions the IV pole at the head of the bed for easy storage.	Yes
DC Motors: Provides quiet and smooth bed operation.	Yes
Frame Options (additional cost)	Hill-Rom Advanta™ 2 Bed
Scale: In-bed scale accurate to +/- 1% of patients weight, even when head is elevated.	Yes
Advanced Bed Exit System: 3 mode, Patient Positioning Monitor allows the ability to change volume and/or tone; ability to alert if system is compromised.	Yes
Powered Transport Device	Yes (IntelliDrive® Powered Transport)
NaviCare® Communication Module	Yes
SafeView™ Alerts: Visually reminds caregivers when bed exit rails, brakes and bed height are in their safe positions.	Yes
Nurse Call, Universal TV, and Lighting Controls	Yes
5" or 6" Urethane Casters with Central Locking Caster System	Yes
Digital Head of Bed Indicator	Yes
Bi-directional EMR Connectivity	Yes
Patient Pendant Control	Optional

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2009 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

161547 rev 1 11/5/09



Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-445-3730
Canada 800-267-2337
www.hill-rom.com