

The Vest® Airway Clearance System

An easy-to-use medical device for children and adults who have excess airway mucus and difficulty clearing it effectively.

- Provides secretion clearance superior to standard chest physiotherapy
- Delivers reliable, technique-independent therapy requiring no special positioning or breathing maneuvers
- Easy to administer and comfortable for the patient
- Chest circumference up to 67"



Art of Care® Large Capacity Chair

- Four unique aesthetic styles
- STEELock™ construction ensures superior structural integrity
- Replaceable components for easy updates
- ConformBack™ option available for extra back comfort
- Fully suspended seating construction supports 750 lbs in 30" or 40" seat width without the use of uncomfortable plywood decking, for superior comfort



Art of Care® Bariatric Recliner

- Three-position seating
- 660 lb weight capacity with 30" seat
- Front and rear locking casters
- Accessory hooks, optional removable arm and tray table



Hill-Rom provides everything you need to care for your bariatric patients with dignity, keep your caregivers and patients safe, and help you manage rising costs.

Our long experience in the care setting means we know how to create products and programs that combine safety with outstanding patient care. We can help you find the right formula for your facility – the right products, the best services, the right mixture of rental vs. capital purchases, the most appropriate ergonomic programs.

Bariatric Services

Environmental Solutions

The weight and body mass of bariatric patients can challenge existing room design. More space, floor mounted toilets, and larger bathing areas are some of the issues that need attention. Hill-Rom's Design Services offers a Bariatric RoomBuilder® program to assist you in designing a room that meets all the unique needs of your bariatric patients.

Ergonomic Solutions

Hill-Rom offers a broad range of consultative and training programs to help protect caregivers who deal with bariatric patients. They include a Detailed Ergonomics Evaluation Process (DEEP), Ergonomics Program Audit, Planning and Policy Development, and Staff Training and Education.

Educational Solutions

Our educational seminar programs provide hospital-based clinicians with practical information for the proactive, comprehensive, clinical management of the severely obese patient. CE programs are available online or through our clinical therapy specialists on a wide range of topics related to the bariatric patient.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2011 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

CTG711 rev 7 5/24/11

USA 800-638-2546
Canada 800-267-2337

www.hill-rom.com



Bariatric Products/Services Guide

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

TotalCare® Bariatric Plus System

One therapy system ... Safety and Efficiency. Mobility and Dignity.

Patient benefits

- Weight-based pressure redistribution
- Available with low air loss with and without Turn Assist and foam options
- Available with foam, low air loss and TotalCare SpO₂RT® surface; integrated air surface provides Advanced Microclimate® technology—the next generation low air loss—for cooler, dryer skin
- Patient helper allows easy repositioning
- Extra wide, 40" surface
- Customizable egress available from foot or side
- Aesthetic design for patient dignity and satisfaction
- Seat deflate for easier patient egress and ingress with low bed height of 21.5" – available with low air loss only

Caregiver benefits

- Turn Assist facilitates easier and more efficient patient handling
- FullChair® position with cradle transition supports the bariatric patient to minimize sliding
- Powered transport for ease of movement throughout the facility
- Accessory outlet at bed site for efficient care
- 30° head-of-bed angle alarm and trending for patient care protocols
- Adjustable length helps keep patients from sliding down in bed
- Weight Limit: 500 lbs

Triflex™ II Bariatric Bed

- Expandable 37-48" width frame
- "Fold and roll" for storage
- Drop down rails
- 3 piece interchangeable foam width mattress
- Pendant and scale included
- Trapeze optional
- Optional Synergy® Air Elite low air loss surface
- Weight Limit: 1,000 lbs



Excel Care® ES Bariatric Bed

Safe, Simple, Dignified.

Patient benefits

- Weight-based pressure redistribution
- Available with low air loss with Turn Assist, low air loss and foam options
- Safely supports patients up to 1,000 lbs
- Robust, easy-to-grip patient-assist siderails
- Aesthetic design for patient dignity and satisfaction
- Expandable width surface offers side air bolster expansion, with both foam and air surface, from 40-50" for greater patient comfort and safety

Caregiver benefits

- Adjustable length helps keep patients from sliding down in bed
- Seat deflate for easier patient egress and ingress – available with low air loss only

Other features

- Low 17" bed height, built-in scale, battery backup, optional patient helper, and optional pressure redistribution, low air loss therapy surface
- Available with Turn Assist for easier and safer patient care

Hill-Rom® Stretchers

Mobility. Comfort. Quality.

- Reinforced frame for all stretchers supports 700 lbs
- Optional, built-in scale shows kilograms and pounds simultaneously
- Available Tempur-Pedic® mattress for patient comfort
- Optional active hand brake to help negotiate ramps and corners without using your body weight and for safe and easy stops
- Available in 26" and 30" widths
- Weight Limit: 700 lbs



AIRPAL® Patient Air Lift

The AIRPAL lift assists caregivers in the transfer of patients from a bed to a stretcher in a variety of care settings.

- Allows patient to be transferred on a cushioned film of air
- Reduces risks associated with patient handling
- Increases efficiency by reducing assistance required for lifting
- Weight Limit: Up to 1,200 lbs



Barton® I-700 Transfer/Chair

Designed to allow a single caregiver to laterally transfer patients safely and easily.

- Safe, dignified lateral transfers up to 700 lbs
- Easy conversion from a flat transfer position to an upright mobile chair
- Adjustable head, foot and arm rest
- Weight Limit: 660 lbs – rental; 700 lbs – capital



Liko® Lift

Innovative lifting solutions for your specialized needs.

- Viking® XL mobile floor lift
- FreeSpan Ultra Twin™ overhead lift
- Very easy to handle
- Ideal for heavy lifting situations, including from the floor
- Weight Limit: 660 lbs (mobile); 1,100 lbs (overhead)



Walker

Our bariatric walker allows patients to ambulate on their own and reduces the need for caregiver assistance.

- Adjustable to accommodate patient height
- Collapsible for easy storage
- Anodized double-lumin tubing for increased strength
- Weight Limit: 800 lbs



Wheelchair

Designed to allow self-mobilization, a key to restoring patient dignity, our bariatric wheelchairs help caregivers get their patients moving.

- Anti-tip mechanism that helps prevent the chair from tipping backwards
- Adjustable, elevating leg supports with calf support
- Full-length removable, height-adjustable arm rests
- Rear push bar provides access for multiple caregivers
- Weight Limit: 800 lbs



Shower/Commode Chair

Addressing both toileting and hygiene needs, our shower/commode chair is sturdy and versatile.

- Adjustable height allows for placement over hospital toilets and easier patient transfer
- Stainless steel construction increases versatility, as chair can be used in the shower
- Rubber stops on all legs to help secure placement
- Weight Limit: 800 lbs

