

# Ampera™ bed

Trusted performance for all care environments



Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**



## Optimal care

Healthcare facilities are faced with the challenge of how to deliver quality care and maintain the dignity of patients/residents. Facilities also require products that are flexible in design and can be combined and configured to provide solutions that assist in the reduction of the risk of falls and work-related injuries, as well as enhance mobility.

At the same time, resource constraints mean that facilities need to ensure maximum value for their healthcare budgets. Institutions require the flexibility of options when considering new capital expenditures.

The Ampera family of beds from Hill-Rom delivers the essential requirements needed in today's care facilities. With advanced features for patient/resident mobilisation and safety, the Ampera range offers a choice when looking at beds that will work with staff to ensure optimal care.



Style



Classic



Select



## Safe and dignified mobilisation

The Ampera features patented **Bed Height Memorisation** (available for Style and Select configurations), that offers the ability to memorise the preferred bed height for an individual. Additionally, the **one-touch egress position** facilitates easier pivoting of the patient/resident to a side-sitting position at the memorised height. This allows for a safer and more comfortable bed egress.

Features such as **Intelligent AutoContour™** and **Sliding Backrest**, are both designed to minimise migration of the patient/resident in the bed and improve comfort, while helping to reduce pressure on the midsection, making breathing easier.

Style and Select configurations of the Ampera also offer **optional retractable egress handles**, which can be raised in one step and secured to provide support for the patient/resident when exiting the bed. When not in use, the handles can be retracted and easily stored.

Also included in the Ampera Select is the Hill-Rom **one-touch EasyChair™** function, which gives the patient/resident the ability to effortlessly reach a sitting position in the bed.



One-touch egress position



One-step egress handles



EasyChair

# Essential patient and resident safety

Falls can lead to loss of autonomy, as well as having other implications for both the care facility and the patients/residents. The Ampera bed is designed to assist caregivers in preventing falls for patients/residents.

This versatile bed can be lowered to a resting position of 35 cm\* to minimise the risk of falling, or raised to a stable working height of 77 cm\* to ensure a safe and ergonomic working position for caregivers. In its lowest position, the bed's proximity to the floor can give fall-prone patients/residents a greater feeling of safety and help to reduce the severity of falls.



Optional full-length wooden siderails



Ampera Select offers an **intelligent night light**, which alerts caregivers when the bed is either in the lowest position (green) or out of the lowest position (amber). It also serves to assist in patient/resident orientation at night.

Ampera Select also incorporates a **brake-off indicator**, and all configurations offer **centralised braking** on all four castors for maximum patient/resident security.



## Caregiver safety and efficiency

### Increased manoeuvrability with an ideal ergonomic design

The ergonomic design of Ampera combines light weight and high stability, reducing the effort required to transport the bed, thereby helping to prevent caregiver injuries.

- Optional **150 mm castors** mean caregivers can manoeuvre the bed with ease
- Steering functions are made easier via the **steering castor** facilitating correct directional movement
- In the raised position the **multi-functional brake** pedal also acts to control the steering castor



## Safe mobilisation for caregivers with Liko™

By combining innovative bed systems and lifting solutions from Liko (a Hill-Rom product brand), along with structured mobility and fall prevention programs, caregivers can mobilise patients and residents safely, easily and more frequently with less risk of caregiver injuries. The Ampera bed, along with Liko mobile and overhead lifts, promote independent mobility and dignified care.

# Designed to optimise choice

The Ampera bed is offered in three distinct configurations, each with a range of corresponding options, combining the necessary care functions in tandem with enhanced styling to improve the comfort and well-being of patients/residents.



## Ampera Classic

- **No siderails** mean greater independence for patients/residents and a more home-like environment
- A **220 kg Safe Working Load** accommodates a greater range of patients/residents
- **Centralised braking on four castors** ensures maximum safety for patients/residents



## Ampera Style

*(includes all features/functions of the Ampera Classic and additional features below)*

- **Patented bed height memorisation** means an accessible, easily pre-programmable height for patient/resident bed egress
- **Egress positioning** contributes to easier and dignified patient/resident mobility
- **3/4 siderails** incorporate a **double lock-out** release, which prevents accidental lowering. They also feature protective covers that can absorb the negative effects of unwanted impacts
- **Optional egress handles** allow for easier and safer bed exit



## Ampera Select

*(includes all features/functions of the Ampera Style and additional features below)*

- **Intelligent night light** increases the feeling of safety by minimising disorientation at night
- **EasyChair function** offers one-touch patient/resident seated positioning
- **Optional bed extension and linen holder** assists caregivers in their daily activities



## A Total Room Solution™ from Hill-Rom

Hill-Rom provides a total integrated room solution to ensure safe, comfortable, and dignified care. Lift systems from Liko work together with Hill-Rom beds to assist patients/residents and caregivers with mobility challenges. Furniture options allow for functional, flexible combinations to complement each bed, including bedside cabinets, overbed tables, wardrobes, chests of drawers, desks, chairs, wall units and shelving.



## Service

Hill-Rom Service—one supplier to meet all your needs

To optimise lifelong performance, we build serviceability into our equipment beginning at the R&D stage. As your dedicated service provider we help you improve productivity and efficiency by optimising the maintenance of your equipment, minimising downtime and repair costs and facilitating compliance with local laws and Quality and Regulatory requirements.

Hill-Rom offers a variety of professional services, performed by certified service technicians, to meet all your maintenance requirements.

Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with healthcare providers in more than 100 countries, across all care settings, by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency and Respiratory Health**. The people, products and programs of Hill-Rom work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**



Hill-Rom is committed to the concept that mobilising people early and supporting their independence, improves life - in the hospital and at home.

The Ampera medical bed is intended for use in residential care facilities and retirement homes for patients with a weight  $\geq 40$  kg, height  $\geq 146$  cm and BMI  $\geq 17$ .

Class I

Manufacturer: Hill-Rom S.A.S. – BP 14 Zi du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: SEN138301-03, 10 July 2017

Not all products/options are available in all countries.  
For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**