

Tri-Flex™ II

Bariatric Bed



Dependable and Economical

- Multiple positions:
 1. Single touch Trendelenburg/ Reverse Trendelenburg functions
 2. Cardiac chair position
- Expandable width frame
- Drop down foot rail
- Independent head and knee positioning
- Scale with easy-to-read display
- "Fold and roll" storage capability
- Foam mattress with two interchangeable 14 cm (5.5") side bolsters
- Optional low air loss
- Optional trapeze
- Impulse Drive System allows for safe and easy control during transport

Tri-Flex II Bariatric Bed

The Tri-Flex II Bariatric Bed is a dependable, economical bed that provides support for both patient and caregiver. Features such as expandable width frame, drop down footrails, battery back up, CPR mechanism and choice of surfaces make this a reliable, basic choice for the care of your bariatric patients.



Bed		Caster and Braking System	
Low position— sleep deck to floor	36 cm (14")	Caster diameter	15 cm (6")
High position— sleep deck to floor	71 cm (28")	Four wheel braking	No
Maximum overall width	104 cm (40.75")	Steer and Lock	Yes
Maximum overall length	239 cm (94")	Scale System	
Patient weight capacity	454 kg (1,000 lb)	Accuracy	1% of patient weight
Bed Angles		Scale weight capacity	454 kg (1,000 lb)
Head section	0-60°	Safety Listing and Power Requirements	
Knee section	0-24°	Safety listing	UL 60601-1
Max Trend	12°	Power requirements	120V AC, 60Hz
Sleep Surface		Standard Features	
Sleep deck width	94-137 cm (37-48")	Expandable width frame and foam surface	
Sleep deck length	218 cm (86")	Pendant	
Foam surface (length x width x height)	218 x 122 x 24 cm (86 x 48 x 9.4") [includes two 14 cm (5.5") interchangeable side bolsters]	Battery backup	
		Optional Features	
		Full frame trapeze	
		In-bed scale	
		Impulse Drive System	

Hill-Rom is a leading global medical technology company with more than 6,500 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom's people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries. For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This product is not available in countries requiring a CE mark. This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments. These medical devices comply with all applicable regulations. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Medical beds manufacturer:

Hill-Rom SAS - B.P. 14 - Z.I du Talhouët - 56330 Pluvigner - France (also European Representative)
Völker GmbH - Wullener Feld 79 - 58454 Witten - Germany
Hill-Rom Holdings, Inc. - 1069 State Route 46 East - Batesville, IN 47006 - USA
Lifts manufacturer: Liko AB - Nedre vägen 100 - 975 92 Luleå - Sweden

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2014 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN126304-01, 21-07-2014

Enhancing outcomes for
patients and their caregivers:

Hill-Rom