

# Bedside cabinet – high

Crafted to provide a welcoming home environment in today's care facilities..



## Features and benefits

- A varied selection of drawer and door configurations, suitable across all budgets
- Soft-close, shock-absorbent drawers that glide effortlessly into place
- Long-lasting, durable construction, using materials that are easy to clean and maintain
- Aptly sized to suit the needs of today's care environments

## Bedside cabinet – high

## Edison Collection - Classic



Optional: Sliding tray



Optional: Accessory bar holder



Optional: Simple base



Optional: Two half doors



Optional: Fixed strip with optional colour contrast



Optional: Colour contrast for fixed shelf

### Technical Specifications model CH700B4

|                                |                              |
|--------------------------------|------------------------------|
| Height                         | 885 mm                       |
| Width                          | 401 mm                       |
| Depth                          | 401 mm                       |
| Weight                         | ~22 – 29.5 kg                |
| <b>Material</b>                |                              |
| Body                           | Melamine                     |
| Top plate                      | High Pressure Laminate (HPL) |
| Edges                          | Aluminum                     |
| <b>Hinge for B01, B02, B03</b> |                              |
| standard                       | on the left side             |
| optional                       | on the right side            |
| <b>Fixed strip</b>             |                              |
| standard                       | for B02, B03, B04            |
| optional                       | for B01                      |

### Castors

|          |   |
|----------|---|
| standard | 4 x 50 mm with brakes                                     |
| optional | simple base<br>(2 x 50 mm with brakes and 2 end fittings) |

### Handles/Colours

see handles and colours sheet

### Additional options

- sliding tray
- accessory bar holder for one or two sides
- two half doors for B01, B02, B03
- interior removable shelf for B01, B03
- contrasting colours for
  - sliding tray and/or fixed shelf
  - fixed strip and/or fixed shelf

Hill-Rom is a leading global medical technology company with more than 7,000 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries. For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Furniture manufacturer:  
Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.  
©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
Doc. No: 5EN150317-01, 18 June 2015

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**