



## Simplify Short Transfers.

HandySheet™ and HandyTube™ are manual assistive devices that facilitate transfers by reducing friction between the patient and the underlying surface. They can be used with both patients of normal weight and very heavy patients. The sliding sheets in the Handy™ series are made of a high-quality, low-friction fabric that slides very easily, even after long-term use.

HandySheet™ is mainly used to facilitate application of a sling in a seated or recumbent position. HandyTube™ is used, for example, to reposition a patient higher up in bed, to turn a patient in bed or for transfers between two beds.

Unique features of the HandySheet™ include pockets that facilitate positioning in a tight-fitting wheelchair, and loops on the longer version of HandySheet™, with which the sheet can be secured to head end of the bed. Both features are unique to the Handy™ series and are registered Liko™ designs.

### Important Features

- **The fabric** is of very high quality, with very low friction even after many washes
- **The long versions** of the HandySheet™ have loops in the upper edge to facilitate the attachment to the bed. The number of caregivers holding the sliding sheet is thereby reduced.
- **The short HandySheet™** has small pockets which facilitate the application in, for instance, a narrow wheelchair.
- **We have developed unique work methods** which facilitate the application of sliding sheets and which are also comfortable for the patient. With our method the patient does not have to be turned in the bed to place the sliding sheet.
- **HandyTube™** is available in three different models, with different application areas.

### HandySheet™

Liko™ HandySheet™ is an excellent aid for applying sling, especially on very heavy patients, pain-sensitive patients or patients with contractures. HandySheet™ is available in two variants; Short and Long, of which the longer model is equipped with handles all the way around. HandySheets™ are delivered in sets of two and are used in pairs.



### HandyTube™

Liko™ HandyTube™ is a slide sheet that is mainly used to reposition a patient higher up in bed, to turn a patient in bed or for transfers between two beds. HandyTube™ is available in three different models, each with different application areas: Short, Long and Wide.



Read the instruction guides for the respective products for more information on application areas and practical use. You can also contact your nearest Liko representative for information and advice. Find our addresses at [www.liko.com](http://www.liko.com). Always remember to work in an ergonomically correct manner. If the task seems heavy for the caregiver, we recommend the use of other lifting and transfer assistive devices in our extensive range of lifts and accessories.

### HandySheet™

Product name	Prod. No.	Size	Measurements (cm)
HandySheet™ Short, without handles	3715015	M	70x70
HandySheet™ Short, without handles	3715017	XL	120x85
HandySheet™ Long, with handles	3716015	M	90x185
HandySheet™ Long, with handles	3716017	XL	140x215

Also available in jumbo packs (10 pcs/package).

### HandyTube™

Product name	Prod. No.	Size	Measurements (cm)
HandyTube™ Short, without handles	3710015	M	70x90
HandyTube™ Long, without handles	3711015	M	90x185
HandyTube™ Wide, with handles	3712015	M	140x140

Also available in jumbo packs (10 pcs/package).

HandySheet™ and HandyTube™ are a medical device class 1 products. HandySheet™ and HandyTube™ have been tested by an accredited testing institute and comply with all requirements for MDD Class I products and the harmonized EN ISO 10535 standard.

Design and Quality by Liko in Sweden.

Liko is quality-certified in accordance with ISO 9001 and its equivalent for the medical device industry, ISO 13485.

Liko is also certified in accordance with environmental standard ISO 14001.



Hill-Rom and Design and Hill-Rom word mark are trademarks or registered trademarks of Hill-Rom Services Inc. Liko and Design and Liko word mark are trademarks or registered trademarks of Liko R&D AB. HandySheet, HandyTube and Handy are trademarks or registered trademarks of Liko R&D AB.



Enhancing Outcomes for Patients and Their Caregivers.™

France \_\_\_\_\_ +33 (0)2 97 50 92 12  
 United Kingdom \_\_\_\_\_ +44 (0)1530 411000  
 Deutschland \_\_\_\_\_ +49 (0)211 16450-0  
 Nederland \_\_\_\_\_ +31 (0)347 32 35 32  
 Italia \_\_\_\_\_ + 39 (0)2 950541  
 Suisse/Schweiz \_\_\_\_\_ +41 (0)21 706 21 30  
 (deutschsprachig) \_\_\_\_\_ +41 (0)21 706 21 38

Österreich \_\_\_\_\_ +43 (0)2243 28550  
 Ireland \_\_\_\_\_ +353 (0)1 413 6005  
 Iberia \_\_\_\_\_ +34 (0)93 6856000  
 Nordic Countries \_\_\_\_\_ +46 (0)20 78 10 30  
 Export \_\_\_\_\_ Contact your local Hill-Rom distributor or contact your Area Manager via website or call \_\_\_\_\_ +1 812 934 8173

[www.hill-rom.com](http://www.hill-rom.com)  
[www.liko.com](http://www.liko.com)

Hill-Rom® reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom® makes is the express written warranty extended on the sale or rental of its products.  
 © 2012 Hill-Rom® SARL, ALL RIGHTS RESERVED.