

Uno™ 102 EE mobile lift

Ideal for use in homecare and care facilities



Features and benefits

Three variable height settings accommodate a variety of different lifting situations

The **FlexLink** (an extra link arm that increases the distance to the lift mast) **facilitates vertical lifting**, as well as offers greater flexibility when connecting the sling to the slingbar

Manual and electric emergency lowering mechanisms ensure safe operation for caregivers

Optional quick release hook system means caregivers can easily change a sling bar or install a scale

The compact design offers a **greater range of access for lifting** throughout the home and care facilities

An extensive assortment of accessories is available including slings, slingbars and scale

Uno 102 EE mobile lift



Flex link for improved patient positioning



Manual emergency lowering



Hand control

Technical Specifications

Maximum Load	175 kg
Weight	42.7 kg
Height	1435 to 2015 mm
Protection class	IPX4
Slingbar options	
standard	Universal 450 (450 mm)
optional	Mini slingbar 220, Universal slingbar 350, Universal TwinBar 670, Crossbars
Base	
height	100 mm
length	1255 mm
width	690 - 1090 mm
Base-width expansion	electric
Castors	
front	75mm double band
rear	75 mm double band with brakes

Material	powder coated steel
Emergency lowering	manual and electric
Lifting speed	30 mm / sec
Maximum lifting height	1790 mm
Lifting interval	1270 mm
Sound level	39 dB(A)
Battery	2 x 12 V 2.9 Ah valve-regulated lead-acid gel-type
Charger	built-in
Hand control	wired
Hand control protection class	IPX 5
Product number	2010004

Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with healthcare providers in more than 100 countries, across all care settings, by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency and Respiratory Health.** The people, products and programs of Hill-Rom work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

The Liko Uno 102 EE is intended mainly for use in Post Acute Care facilities like nursing homes and other care homes for most common lifting situations.

Class: I

Manufacturer: Liko AB, Nedre vägen 100, SE-975 92 Luleå, Sweden

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN183327-01, 21 April 2017

Enhancing outcomes for patients and their caregivers:

Hill-Rom