

# Wood armchair – SB404A

Seating options designed for today's care environments



SB404A

## Features and benefits

- An extensive selection of fabric and colour configurations, ideal for a range of facilities
- Long-lasting, durable construction, using materials that are both easy to clean and maintain
- Aptly sized to suit a variety of needs across most healthcare settings

## Technical Specifications model SB404A

<b>Height</b>	860 mm
<b>Seat height</b>	490 mm
<b>Width</b>	580 mm
<b>Depth</b>	590 mm
<b>Weight</b>	~8 kg
<b>Maximum load*</b>	135 kg
<b>Material</b>	
fabric	coated polycotton blend
foam	polyurethane
frame/panel frame	plywood
seat frame	solid wood
<b>Foam depth</b>	
seat	45 mm
back - front	40 mm
back - back	6 mm
<b>Fabrics/colours</b>	
see fabrics and colours sheet	
<b>Additional options</b>	
2 castors on front legs	

\*as per manufacturer's recommendation

Hill-Rom is a leading global medical technology company with more than 10,000 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health**. Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products. ©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. Doc. No: 5EN150339-01, 09 October 2015