

Metal recliner – SM608A

Seating options designed for today's care environment



Features and benefits

- Choose from an extensive selection of fabrics, ideal for a range of facilities
- Dual rear castors allow for greater manoeuvrability
- Height-adjustable armrests and retractable headrest enhance comfort and safety for patients
- Body contouring seat, adjustable back and retractable tilting headrest offer increases patient/resident comfort
- Removable back and seat cushions allow for easier cleaning of the chair



Head support



Comfort seating pad



Overchair table

Technical Specifications **model SM608A**

Height	1,265 mm
Seat height	430 mm
Width	680 mm
Depth	690 mm
Weight	25 kg
Maximum load	130 kg
Material	
fabric	coated polycotton blend
foam/armrests	polyurethane
frame/panel frame	epoxy-coated light grey metal
Foam depth	
seat	60 mm
back	70 mm

Fabrics

see fabrics sheet

Accessories

- head support
- comfort seating pad
- overchair table

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Furniture manufacturer:
Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN150351-01, 10 February 2016