

# Ampera™ bed – Classic

An affordable care option



## Features and benefits

The **lower height of the bed** assists in minimising the risk of fall-related injuries

A **220 kg Safe Working Load** accommodates a greater range of patients/residents

**Central braking and steering functions** ensure enhanced security for patients/residents

**Line-Of-Site™** backrest indicator ensures accurate patient/resident positioning information for caregivers

**Intelligent AutoContour™ and Sliding Backrest** minimise patient/resident migration in bed, reducing abdominal compression during bed articulation

**Caregiver control with total lock-out** means enhanced operational safety

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**

## Ampera bed - Classic



Pendant control



Bilateral mechanical CPR



Central braking



Optional: Full-length wooden siderails

### Technical Specifications (for standard configuration)

|   |                    |
|---|--------------------|
| Safe Working Load                               | 220 kg             |
| Bed weight                                      | 137 kg             |
| <b>Height adjustment</b> (dependent on castors) |                    |
| standard  | 35 to 73 cm        |
| optional  | 39 to 77 cm        |
| Protection class                                | IPX4               |
| <b>Bed angles</b>                               |                    |
| Head section                                    | 0 to 65°           |
| Thigh section                                   | 0 to 28°           |
| Foot section                                    | -3 to -22°         |
| Sleep deck width                                | 90 cm              |
| Sleep deck length                               | 200 cm             |
| Sleep deck surface                              | metal              |
| <b>Siderails</b>                                |                    |
| standard  | none               |
| optional  | 3/4 length metal   |
| optional  | full-length wooden |

### Castors

|          |                           |
|----------|---------------------------|
| standard | 4 x 100 mm<br>single band |
| optional | 4 x 125 mm<br>single band |
| optional | 4 x 150 mm<br>single band |

### Steering castor

|          |             |
|----------|-------------|
| standard | at head end |
| optional | at foot end |

### Central brake

at foot end

### Lifting/IV pole holder

inside double holder  
at head end

### Bilateral mechanical CPR (Cardiopulmonary Resuscitator)

### Controls

see controls sheet

### Designs/Colours

see designs and colours sheet

Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with healthcare providers in more than 100 countries, across all care settings, by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency and Respiratory Health.** The people, products and programs of Hill-Rom work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

The Ampera medical bed is intended for use in residential care facilities and retirement homes for patients with a weight  $\geq 40$  kg, height  $\geq 146$  cm and BMI  $\geq 17$ .  
Class I

Manufacturer: Hill-Rom S.A.S. – BP 14 Zi du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
Doc. No: 5EN138302-04, 10 July 2017

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**