

# OctoStretch™

Maximise horizontal lifting space



With a lifting capacity of 200 kg, OctoStretch is ideal for patient transfer in a variety of healthcare settings, including surgical, radiology and Intensive Care Unit (ICU) departments. OctoStretch can be used in combination with all Liko overhead lifts, as well as a selection of mobile lifts.

## Features and benefits

**Optional Stretch Leveller™** allows for quick adjustment of the centre of gravity during patient lifting

Can be used in conjunction with either a **full lift sheet** or with **four lift straps**, meaning easier access when performing patient care tasks

At just 7.2 kg, the light weight allows for **easy operation** of the lift stretcher

**Compact design** for easier transportation and storage

**Optional quick-release hook system**, which means caregivers can easily and rapidly facilitate transfers between many lift systems

# OctoStretch



Optional Stretch Leveller

## Technical Specifications

Maximum Load	200 kg
Weight	7.2 kg
Size	
height	227 mm
length	1263 mm
width	617 mm
Material	powder-coated steel
Product number	3156055
with stretch leveller	3156056

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

OctoStretch accessory is intended for lifting and transferring patients in a horizontal position in cases where a particularly flat underlying lifting surface is not required. Common areas of use are in surgical wards, however it is also used in intensive care wards and X-ray departments.

Class I

Manufacturer: Liko AB, Nedre vägen 100, SE-975 92 Luleå, Sweden

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
Doc. No: 5EN183316-01, 07 February 2017

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**