

# Hill-Rom® Surgical Stretcher



## Go Beyond Mobility™

With its 700 lb weight capacity, the Hill-Rom Surgical Stretcher provides flexibility to comfortably accommodate all patient conditions as well as your needs. Its simple-to-use articulating head section gives surgeons optimum positioning control to perform surgical procedures.

A rush to the operating room is not just about running down hallways and turning corners. It's about the minute details that contribute to a safe and comfortable experience for both you and your patient.

Our stretchers **provide more than mobility** for your patients:

- A safe and comfortable patient experience
- Simple and safe maneuvering
- Durability to withstand your demands

## Key benefits at a glance:

### A Safe and Comfortable Patient Experience

- Help prevent pressure injuries with the pressure redistribution ComfortLine® Mattress
- Provide extra room for recovery with PACU extenders
- Protect your patient's hands with bumper protected siderails
- Provide a gentle patient ride — with Steering Plus™ system that absorbs the shock of driving over bumps and thresholds
- Allow for safe patient ingress and egress with low height

### Simple and Safe Maneuvering

- Optimum positioning control through articulating head section
- 360° Maneuvering with Steering Plus™ system
- Zero Transfer Gap — with OneStep® Tuck-Away Siderails

### Durability to Withstand Your Demands

- 700 lb weight capacity
- Compatibility with most lifts and overbed tables
- Easy to clean smooth deck and siderail design
- Power washable

# Go Beyond Mobility™

A safe and comfortable experience for both you and your patient starts on a Hill-Rom® stretcher.



Articulating Head Section



Steering Plus™ System



Temporal Wrist Rest



Superior Wrist Rest

A SAFE & COMFORTABLE PATIENT EXPERIENCE	STANDARD FEATURES	OPTIONAL FEATURES
26" width	X	
Manual Knee Gatch	X	
Bumper Protected Siderails	X	
Trendelenburg/Reverse Trendelenburg	X	
Six Transport Strap Locations	X	
PRESSURE REDISTRIBUTION MATTRESS OPTIONS		
4" Surgical Surface with Flat/Concave Head Cushion Set		X
4" Surgical Surface with Flat Head Cushion		X
5" AccuMax Quantum™ mattress (Elastimax or Hypolex)		X
SIMPLE & SAFE MANEUVERING		
Steering Plus™ System	X	
OneStep® Tuck-Away Siderails	X	
BackSaver® Fowler feature		X
Brake/Steer – all four sides		X
Dual-Sided hi/lo Hydraulic Control Pedals	X	
Ergonomic Push Handles	X	
Ergonomic Push Handles with Integrated IV Pole Transport Holder		X
3/4" Standard IV Pole Sockets	X	
Removable IV Pole – 2 Stage		X
Permanent IV Pole – 3 Stage (40 lb weight capacity)		X
Radiolucent Deck	X	
Upright Chest X-Ray Cassette Holder	X	
DURABILITY TO WITHSTAND YOUR DEMANDS		
700 lb weight capacity	X	
Full Perimeter Bumping	X	
Smooth Deck Design	X	
Decals Options	X	
Bumper Color Options – choice of neutral, blue, red, teal or purple (green bumpers only for IntelliDrive® models)	X	

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
198177 rev 1 29-JUN-2016 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730  
Canada 800-267-2337

[www.hill-rom.com](http://www.hill-rom.com)

Technical Specifications	
Maximum weight limit	700 lbs (317.5 kg)
Low position (floor to top of sleep deck)*	21.5" (54.6 cm)
High position (floor to top of sleep deck)	33.25" (84.5 cm)
Overall length	92" (233.7 cm)
Overall width (siderails up)	32" (81.3 cm)
Overall width (siderails stored)	29.25" (74.3 cm)
Siderail length	48" (121.9 cm)
Siderail height above sleep deck	14.5" (36.8 cm)
Comfortline® mattress size	30" (76.2 cm) x 75" (190.5 cm)
Maximum head elevation	90°
Maximum Trend/Reverse Trend	18°
Floor to base clearance	3.5" (8.9 cm)
Caster size	8" (20.3 cm)
Warranty	
Frame & Welds for life of product	
Three-year parts	
Two-year mattress	
One-year labor	
Accessory Options — Please confirm accessory availability for model of stretcher ordered	
Accessory IV pole	Liquid oxygen tank holder
Ankle stirrups	Oxygen tank holder
Armboard	Paper roll dispenser
Chart holder	Patient tray
Convertible footboard	Restraint straps
Detachable footboard	Siderail covers
Foot extender	Stretcher pillow
Infusion support system	Transport straps
IV transporter	Utility tray
Lateral X-ray cassette holder	

\*Dependent on configuration

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**